

New Product Introduction



HOLD MY CPAP HOSE™ CPAP Tubing Hose Holder



Item	HOLD MY CPAP HOSE™ CPAP Tubing Hose Holder
Item #	BJ240182
Features	<ul style="list-style-type: none"> • Prevents 'fighting', 'tug-of-war' and 'neck choking' of the CPAP tubing while sleeping • Helps to prevent CPAP mask leaks as you toss and turn • Top tubing and adjustable arm loop swivels and moves as you move - allowing the user to move around at night without getting tangled in the CPAP hose • Elastic internal cord stretches as the user moves (figure A) • Base slides between box spring and mattress - installs in seconds, no tools required! (figure B) • Folds flat for easy travel and portability • Includes a travel carry pouch

