


## New Product Introduction



# TAKE THE PRESSURE OFF™ Deluxe Adjustable Pressure Pump and Pad with Flaps



Item	TAKE THE PRESSURE OFF™ Deluxe Adjustable Pump & Pad with Flaps	
Item #	BJ145720	
Features	<p><b>Bubble Pad with Flaps</b></p> <ul style="list-style-type: none"> <li>• For both hospital and home care use</li> <li>• A-B lines alternate inflating 96 air bubbles for superior bedsore treatment (A)</li> <li>• Designed to be placed on top of an existing mattress</li> <li>• 15.7" end flaps wrap underneath mattress for increased stability (B)</li> <li>• Durable, heavy gauge, premium medical grade PVC</li> <li>• Inflated pad size: 78.7" L x 35.4" W x 2.8" H</li> <li>• Tubing and repair kit included (C)</li> </ul>	<p><b>Alternating Pressure Pump</b></p> <ul style="list-style-type: none"> <li>• Provides variable pressure setting for optimum comfort and pressure</li> <li>• Produces 6 LPM of consistent air flow and pressure</li> <li>• Preset 12 minute alternating cycle</li> <li>• Pump features built in brackets in order to easily attach to the end of the end</li> <li>• Quiet pump alternately inflates and deflates 96 individual air cells in the bubble pad</li> <li>• Lighted ON/OFF switch</li> </ul>
Barcode #	652508380308	