

## New Product Introduction



# BUCKLE UP™ Wheelchair Seat Belt with Auto-Type Closure



Item	BUCKLE UP™ Wheelchair Seat Belt with Auto-Type Closure
Item #	BJ270100
Features	<ul style="list-style-type: none"> <li>• Helps prevent falling or sliding out of a wheelchair - great for elderly, epileptic or those patients with loss of muscle function</li> <li>• Easily attaches to most wheelchairs or transport chairs</li> <li>• Buckle is simple to fasten, and releases easily with the push of a button just like a vehicle or airplane seat belt</li> <li>• Strong and durable</li> <li>• Belt measures 1.9" wide with an adjustable length from 35" to 53"</li> </ul>

