

New Product Introduction



TAKE A BREATH™ Mr. Blue Bear Pediatric Compressor Nebulizer



Children Love
Mr. Blue Bear!



Includes 2 child masks, medication cup,
tubing, t-piece, mouthpiece, and 5 filters.



"Hibernation" Carry Bag

Item	TAKE A BREATH™ Mr. Blue Bear Pediatric Compressor Nebulizer
Item #	BJ240400
Features	<ul style="list-style-type: none"> • The child-friendly design reduces anxiety and makes medication time fun and easy • Converts prescribed medication into an aerosol mist of easily inhaled microscopic droplets • Durable, high-power compressor system can dispense all types of liquid medication quickly and effectively • Delivers 39 psi at maximum pressure and 16 psi operating pressure (liter flow range 6-8 LPM) • Operates quietly below 40 dBA • Lightweight and easy to carry with a convenient built-in carrying handle • Comes with two masks, suitable for a wide range of both infant and child users
Warranty	5 year warranty
Barcode #	652508380407



Retail box