

New Product Introduction



WALK WITH ME™ Aluminum Adjustable Crutches



Retail packaging



The handgrips are easy to move and adjust



Push-pin adjustment makes it easy to adjust the crutch height in 1" increments

| | |
|----------|--|
| Item | WALK WITH ME™ Aluminum Adjustable Crutches |
| Item # | BJ210184 (Youth Size) BJ210180 (Adult Size) BJ210188 (Tall Size) |
| Features | <ul style="list-style-type: none"> • Aluminum crutches are durable, lightweight and provide stability • Crutch design provides additional support to the weight-bearing area • Push-pin adjustment makes it easy to adjust the crutch height in 1" increments • Comfortable, durable accessories are latex free • The large-size crutch tips provide excellent traction • The handgrips are easy to move and adjust • Weight capacity: 350 lbs. • BJ210184 for user height 4'6" - 5'2" • BJ210180 for user height 5'2" - 5'10" • BJ210188 for user height 5'10" - 6'6" |



Crutch design provides additional support to the weight-bearing area



The large-size crutch tips provide excellent traction