

New Product Introduction



BATHROOM PERFECT™

5" Raised Toilet Seat with Lock and Removable Padded Arms



Case of 2



Lock keeps the raised toilet seat securely in place

Item	BATHROOM PERFECT™ 5" Raised Toilet Seat with Lock and Removable Padded Arms, Case/2
Item #	BJ110210
Features	<ul style="list-style-type: none"> • Fits most standard and elongated toilets • Adds 5" of height to your toilet • Lock keeps the raised toilet seat securely in place • Padded arms provide balance assistance • Rubber pads help prevent the raised toilet seat from slipping off the toilet • Tool-free installation • Weight capacity: 300 lbs./136 kg. • Width between arms: 17" • Dimensions: 20" wide x 16" deep • Height of arms above the raised toilet seat: 6" high
Item #	652508380889

Bathroom Perfect™ 5" Raised Toilet Seat with Lock and Removable Padded Arms

Item # BJ110210

- Fits most standard and elongated toilets
- Adds 5" of height to your toilet
- Lock keeps the raised toilet seat securely in place
- Padded arms provide balance assistance
- Rubber pads help prevent the raised toilet seat from slipping off the toilet
- Tool-free installation



Retail box