

New Product Introduction



Offset Handle Cane with Deluxe 360° Pivoting Cane Tip



4 Beautiful Colors to Choose From

| | |
|----------|--|
| Item | Offset Handle Cane with Deluxe 360° Pivoting Cane Tip |
| Item # | BJ210190B (Sapphire Blue) BJ210190G (Emerald Green) BJ210190PR (Amethyst Purple) BJ210190R (Ruby Red) |
| Features | <ul style="list-style-type: none"> • Manufactured with sturdy, aluminum tubing • Easy-to-use, one-button height adjustment with locking ring • Handle height adjusts from 30" to 39" • The flexible cane tip pivots to maintain contact with the ground when the cane is at an angle, providing excellent balance • 3 point cane tip allows the cane to stand on its own • Ergonomically designed handle with soft foam rubber grip provides comfort and security • Integrated wrist strap • Weight Capacity: 300 lbs. |



360° Pivoting Cane Tip



ALSO AVAILABLE: BJ210101 - Pivoting Cane Tip ONLY