



New Product Introduction



PERFECT MEASURE™ Deluxe Blood Pressure Kit with 2 Cuffs (Standard and Large) and Accessories



Item	PERFECT  MEASURE™ Deluxe Blood Pressure Kit
Item #	BJ120106
Features	<ul style="list-style-type: none"> • Kit Includes: (1) blood pressure monitor, (1) standard adult arm cuff - 8.6" - 12.6", (1) large adult arm cuff - 8.6" - 16.5", (4) 'AAA' batteries, (1) USB type-C charging cable, (1) user manual • Accurate results detected by oscillography • Fully automatic one-touch operation • Bright, easy-to read LCD screen • Four color backlit display • 2 user mode • 198 total memory recall (99 memories for each user) • Blood pressure classification indicator • Advanced averaging of 3 most recent readings • Irregular heartbeat detector • Automatically powers off after one minute of non-use to prolong battery life • Operates on four 'AAA' batteries or USB type-C cable (included) • 2 year warranty
Barcode #	652508381213

