


New Product Introduction

**In-Stock
for Immediate
Delivery!**



PERFECT ♥ MEASURE™

Elite Fully Automatic Blood Pressure Monitor with Universal Adult Cuff



Item	PERFECT ♥ MEASURE™ Elite Fully Automatic Blood Pressure Monitor with Universal Adult Cuff	
Item #	BJ120050	
Features	- Includes: (1) blood pressure monitor, (1) adult arm cuff, (4) 'AAA' batteries, (1) USB type-C charging cable - Fuzzy logic algorithm - Fully automatic operation - Accurate measurement - Two user mode 180 total memories (90 memories for each user) - WHO blood pressure classification indicator - Irregular heartbeat indicator - Advanced averaging of the 3 most recent readings - Cuff check indicator - Movement error indicator - Can speak your BP and pulse rate results (optional) - Low battery indicator - Automatic power-off to conserve battery life - Bright, backlit LCD display screen - Easy-to-read display shows systolic, diastolic and pulse rate readings - Universal adult cuff fits arms 8.7" to 16.5" in circumference 2 power source options: 4 'AAA' batteries or USB type-C charging cable (included) 3 year warranty	
Barcode #	652508381572	

